



Custom
Solutions
Group

Custom Solutions Group LLC
1419 Avenue A, Katy, TX 77493
Office: 281.574.9999
E-Mail: sales@csg-labs.com



Agilent Technologies

Channel Partner

Description:

The system features ASTM Method D4815, "Standard Test Method for Determination of MTBE, ETBE, TAME, DIPE, and C1 to C4 Alcohols in Gasoline by Gas Chromatography."

Peak Identification

1. Methanol
2. Ethanol
3. IPA
4. t-Butanol
5. n-Propanol
6. MTBE
7. i-Butanol
8. ETBE
9. TAA
10. DME
11. n-Butanol
12. Benzene
13. TAME




Figure 1, Chromatogram for ASTM D4815

Features and Benefits

- Separate temperature zone for TCEP micro-packed column
- Simple design, rugged construction
- Best serviceability
- Decades of experience
- Made-in-the-U.S.A.

